VisualAITM for Sample-Agnostic Flow Imaging Microscopy Particle Analysis





Abstract

- well as those from mixed samples

Materials and Methods

Protein aggregate prep:

- Prepared ~1 mg/mL solutions of BSA, ovalbumin, IVIg, or NIST mAb in PBS
- Stressed 5 mL aliquots of each protein formulation via:
- Freeze-thaw: 4 cycles x (freeze @ -20°C, thaw @ 25°C)
- Shake: Plate rocker at max speed for 4 hours
- Silicone oil emulsion prep: Ο
 - Blended 5% silicone oil in PBS mixtures for 20 s
- Mixture preparation: Ο
 - Prepared 1.5 mL samples containing 25%, 50%, and 75% BSA aggregates and the remainder silicone oil droplets
- Flow Imaging Microscopy (FIM)
 - Analyzed each protein sample + one emulsion on a different FlowCam 8100
 - Varied flow rates and focus settings to simulate use cases

Results

- Pure Sample Performance (Figure 1):
 - >90% classification accuracy on particles >3 µm in size **regardless of the protein** and instrument used
- Similar accuracy on an out-of-focus unit
- Mixture Performance (Figure 2):
 - Reported protein aggregate content is consistent and linear with actual content

Austin L. Daniels, Emma M. O'Riley, Austin B. Coco, Sigrid C. Kuebler Yokogawa Fluid Imaging Technologies

• Al approaches are effective for identifying particles in flow imaging microscopy (FIM) data • Al-based image analysis tools rely on large image libraries captured on the FIM instrument(s) and protein formulations to be analyzed. These data requirements can make AI tools difficult to develop • VisualAI is a pre-trained AI utility that classifies images of protein aggregate and silicone oil droplet particles larger than 3 μm in any protein formulation captured on any FlowCam 8100/FlowCam LO • We demonstrate VisualAI's performance on classifying FIM images of protein aggregates and silicone oil droplets as



Figure 1: Sample FIM images of protein aggregates and silicone oil droplets captured on three FlowCam 8100 units. The aggregated protein is denoted below the images. Also shown is VisualAl's classification performance (bottom row) as a confusion matrix.



	FlowCam 8100 #3							
	Out of focus							
	3.26 3.26 3.26 3.30 3.32 3.32 3.70 3.70 3.73 3.73 3.73 3.73 3 4.22 4.22 4.24 4.32 4.34 5.00 5.05 5.10 5.16 5.17							
	3.14 3.1 3.42 3.4 3.61 3.6 4.08 4.	.6 3.18 2 3.44 5 3.66 .10 4.10	3.18 3.1 3.44 3.4 3.68 3.6 4.11 4	19 44 58 • 17				
r		Protein	Silicone Oil	Other				
)	NIST mAb	95%	3%	2%				
	Silicone Oil	8%	91%	1%				



VisualAl Classification View							
File	Edit	Class	sificatio	on S	ort	Show	Pre
	Zo	om In	Zoom	Out	Zoon	n Home	Sł
Prote	ein-Sili	cone O	il Sampl	e			
Prot	ein	Silicon	e Oil	Othe	r		
3.2	1 3	.32	3.36	3.39	9 3	0 .43	3.43
4.8	 8	5.28	5.28	5	.45	5.45	
6.9	1	7.81	8.07	7	8.12	8.	18

Cla Apj	ssification proach	VisualAl	Prev. Al- Based Methods	Property- Based Methods	
		Use pre- trained Al tools to classify images	Develop and apply AI tools to classify images	Use particle properties to classify images	
rmance	Analysis Time	Minutes	Minutes	Hours	
Perfo	Accuracy	High	Very High	Low	
stness	Instrument	High	Low	Moderate	
Robus	Sample	High	Moderate	Moderate	
of-use	Images Needed	None	10,000s	100s	
Ease-c	Software Integration	Yes	No	Yes	



Scan this QR code to learn more about this study, VisualAI, and FlowCam

