



## VisualSpreadsheet Satellite Version

Analyze your FlowCam data on up to 3 external computers



### A satellite license purchase includes:

- The latest version of VisualSpreadsheet for your FlowCam
- Licenses for up to 3 additional computers

### Compatible computer specifications:

- Windows 10 or 11 (64-bit) operating system
- Intel i5 or AMD Ryzen 5 processor or higher
- 512 GB SSD hard drive or higher
- 16 GB of memory or higher

### With a satellite license, you can:

- Work with FlowCam data in its native, VisualSpreadsheet environment from the comfort of your office or home computer
- Analyze data externally, freeing up your FlowCam instrument for new data acquisition
- Combine data from multiple FlowCam instruments into a single database
- Take advantage of VisualAI™ to automatically differentiate protein aggregates from silicone oil droplets in biopharmaceutical applications
- Import data collected using older FlowCam models into the newest version of VisualSpreadsheet

### Contact us to learn more!

